



Please Scan the QR code
to re-view the User Manual

Phone: 0086-4000078388

Service WeChat: migoosmart

Working Hours: 9:00-18:00, Monday to Saturday
(China time, Day off in Holidays and
weekend)

Company: Shenzhen Migoosmart Co., Ltd.

Address: Rm.807, Building U6, Junxiang U8 Intelligent
Manufacturing Industrial Park, Gushu Community,
Xixiang Sub-District, Bao'an District,
Shenzhen, China

E-mail: service@migoosmart.com

Post Code: 518000




WeChat
Official Account



Service WeChat

※ The final interpretation right belongs to Shenzhen Migoo Smart Co., Ltd.





inFace RF Beauty Instrument

U s e r M a n u a l

Welcome to use inFace

inFace RF (Radio Frequency) Beauty Instrument is independently developed and manufactured by Shenzhen inFace Technology Co., Ltd., which is affiliated to Shenzhen Migoo Smart Co., Ltd. Since its inception, the company has always been based on the practicality, functionality and safety of the products. With the aim of creating scientific and smart products, we are striving to explore and innovate.

Function Introduction

inFace RF Beauty Instrument is a safe, effective and convenient professional beauty instrument that delivers real, visible and instant effects. As a household beauty instrument, the instrument presents its radio frequency beauty technology that is comparable to a professional beauty salon. The inFace RF Beauty Instrument focuses on using the 6-pole low-power radio frequency to penetrate the dermis layer, position tissue heating and stimulate the dermal activity to achieve followings: 1. Acting on the dermis layer of the skin to shallow the wrinkles; 2. Promote the regeneration of subcutaneous collagen and generate new collagen to reduce wrinkles, keep skin firm and supple, and delay skin aging effectively.



* Notice: Please read the manual carefully before using the product.

Safety Instructions

All warnings and safety precautions must be observed before using any electronic equipment, all precautions should be read to ensure safety. Failure to follow the instructions in this manual may cause burns or damage to your skin.

Warnings

Important!

1. Do not use the instrument for eyelids, lips, ears, etc.
2. It is forbidden to use the instrument for purposes other than the designated functional care.
3. The instrument is only charged with original charging cable.
4. If the instrument has been disassembled by yourself, do not continue to use it.
5. inFace RF beauty instrument particular gel is for external use only.
6. Keep out of reach of children.
7. Keep your hands completely dry when operating the instrument.
8. inFace RF beauty instrument and the particular gel must be used under the instructions.

Unsuitable for the following conditions:

1. Under 18 years old
2. People with heart disease, coronary heart disease, thyroid disease, nervous system disease, high blood pressure, diabetes, glaucoma, sensitive to current;
3. People currently or have had cancer, especially skin cancer or precancerous;

4. People who are equipped with electronic devices such as hearing aids, pacemakers or internal defibrillators;
5. People in severe concomitant conditions such as heart disease;
6. Pregnant woman, lactating or menstruating women;
7. People have had a tumor, abnormal wound healing;
8. People under hormone therapy;
9. People with impaired immune system caused by immunosuppressive diseases, such as AIDS or are taking immunosuppressive drugs;
10. Please consult your doctor after any surgery;
11. There is any heat-stimulated disease in the care site, such as recurrent simple herpes virus;
12. There is any skin disease such as sores, sputum, rash, etc. in the care site;
13. There are serious lesions, scar hyperplasia, inflammations and other local lesions in the care site;
14. There is severe sunburn in the care site;
15. There is laser rejuvenation, chemical rejuvenation or wound healing in the care site within three months;
16. There was permanent filler in the care site, such as buried metal, plastic, silica gel, etc.;
17. There was temporary filler in the care site, such as hyaluronic acid and botulinum. Please consult your doctor first;

Notice: If you have any doubts about any potential side effects, please consult customer service or a doctor before using the inFace RF beauty instrument.



Cautions

Important!

1. Do not use the instrument for thyroid or lymphatic areas;
2. Do not use the instrument for throat, eyeballs, ears, upper eyelids, temples;
3. Operate carefully when use the instrument for the areas of thin skin and with blood vessels distributed, such as the corners of the eyes, forehead, neck, mouth, and tibia;
4. people with sensitive skin is recommended to use the low level;
5. Within one hour before using the instrument, do not use flammable liquid skin care products such as alcohol, acetone, perfume, disinfectant, etc.;
6. If you feel any discomfort during use, stop immediately and consult a doctor;
7. Do not clean the instrument with wet cloth or liquid directly, wipe it gently with a soft, dry, lint-free cotton cloth or alcohol pad;
8. When operating the instrument, the pole should be in close contact with the skin;
9. Different gears can be selected, higher position delivers higher power and better effect; select appropriate position in accordance with your skin temperature tolerance;
10. Gel should be used in the care site. Using the instrument on dry skin directly will result skin burns. Continuously add gel during the operation. If the gel on skin is not enough, the effect will reduce;
11. It is recommended to wash your face with cold water after the operation;
12. Before operating, take off your contact lenses, metal ornaments and dentures;
13. Do not store or use the instrument in places with strong magnetic fields and electric waves, in places with high humidity and dust, such as bathrooms and kitchens, or in high-temperature places. Storing the instrument under extreme temperatures and conditions may

cause damage to the material of the machine body and reduce its service life;


14. Do not strike, knock or slam the instrument;

15. After using, please clean the instrument in time. When cleaning the instrument, turn off the power first. Alcohol cotton is recommended to clean the instrument, do not use organic solvents such as benzene or thinners to clean the instrument.

Accessories and Product Appearance

The following items should be found in the package:

- ① Machine Body
- ② Charging cable
- ③ User's Manual and Warranty Card




(If any above accessories is missing, please contact our after-sale customer service center)

Product Specifications


Product Name	inFace RF Beauty Instrument
Product Model	MS6000
Product Size	154mm X 63mm X 39mm
Net Weight	About 125g
Battery Capacity	1500mAh
Working Voltage	DC 3.7V
Working Current	DC3.7V_1000 mA
Rated Power	8W
Standby Current	<3uA
Working Frequency	1M Hz
Charging Mode	USB DC 5V
Charging Current	600mA
Charging Time	About 2.5h
Continuous Working Duration	<60 Minutes
Working Temperature Range	0°C~40°C
Material/ingredient	Titanium Head, ABS, PC, Silicone

Operation Method

1. Long press the button for 1.5 seconds to turn on with a beep sound.
2. Short press to switch three strength levels, high-medium-low, with beeps and white LED indicator
3. In the on state, long press the button for 1.5 seconds to turn off, with beeps.
4. With intelligent timing function, the instrument can automatically shut down after 10 minutes when it is on.
5. When the instrument is on and the titanium heads don't touch skin for 10 seconds, the three LEDs will blink for 30 seconds and then the instrument will automatically turn off.
6. With automatic memory function, the instrument can automatically store the last usage level as the default for the next startup.



Long press the button



Three level adjustable

Tips

1. Clean your face before using the instrument.
2. Gel should be used with the RF Beauty Instrument. Before use, the care site should be covered with gel. Do not use the instrument on dry skin directly.
3. Keep the pole in close contact with your skin and slide it at a constant speed, at this time, you will feel the heat output to the depth of the skin. The slower the sliding speed, the stronger the heat feeling. You can adjust the sliding speed according to your feeling. During the operation, keep six poles contact with your skin at the same time. Follow the principle of going from bottom to top, from inside to outside. Slide the pole head in circles on the face, repeat 5 times.
4. Three different levels can be selected, switch it in accordance with your skin adaptability. Higher power leads to better effect.
5. After the care, wash the remaining gel on your face thoroughly.
6. After the treatment of high temperature, the bottom layer of skin needs to be soothed and calmed. It is recommended to use a moisturizing mask and pay attention to daily sun protection.

Using Method



Moving from the upper edge of the collarbone to the lower edge of the mandible (avoiding the throat)



Circling from the mandible to the ear



Circling from the corner of the mouth to the front of the ear



Circling from the decrease to the cheekbones



Circling from the side of the nose to the cheekbones (avoid the temple)



Sliding from the inner corner of the eye to the outer corner of the eye, and raising the end of the eye, the end of the eyebrow (avoiding the eyeball and upper eyelid)



Raising the corner of the eye from the lower eyelid to the eyebrow, hairline



Pulling from the eyebrows (eyebrow heads, eyebrow peaks, eyebrow tails) to the hairline

Using Frequency


1. Light mature skin is recommended to use once every 3-4 days.
2. For mature skin such as rough skin, dry and aging, it is recommended to use once every 1-2 days.
3. For daily skin care with good skin quality, it is recommended to use once every 1-2 weeks.

Cleaning Method

Long press the button for 1.5 seconds to turn off and clean the instrument with alcohol cotton or soft towel.

Charging Method

1. Connect the supplied USB cable with the standard mobile phone adapter (DC 5V) and the micro USB port of the device to charge.
2. One beep sound indicates going into charging state, and the white LED turns on and flows.
3. When the battery is full, three white LEDs are always on, and the white LEDs goes out after the charging cable is unplugged.
4. The button operation is invalid during charging.



Charging Interface

Preservation Method

- 1.If it is not used for a long time, please wipe up the product and store it in a dry and ventilated place.
- 2.In case of danger, do not expose children to the product.

Common Questions

If you find any issues during use, please check the abnormality according to the following table. If you cannot solve the problems, contact the after-sales service.

Phenomenon	Solutions
1 Feel stinging after using the instrument	Before operating the instrument, gel must be used as a medium for the instrument to contact your skin. Without the gel, you will feel obvious stinging and it will cause negative effects on your skin. There are three strength levels, the low level is recommended for the first time. For the following operation, you can switch the levels according to your skin condition.
2 The skin turns red and glows after using the instrument	The RF instrument generates a magnetic field to the skin with a high-frequency RF pulse current, penetrates into the dermis layer, generates heat, and promotes collagen reorganization and regeneration when the thermal energy of the dermal tissue reaches 45 ° C - 50 ° C. It is normal for the skin to become red and slightly swollen after using the instrument, and it fades within 30-60 minutes.



3 Automatically shut down after 10 minutes of use	The instrument comes with 10 minutes of intelligent timing protection function. If there is no operation for 10 minutes, the instrument will automatically enter the protection mode and shut down. You can turn on the instrument again to keep using.
4 Failure to turn on	Long press the button for 1.5 seconds to turn on or turn off. If there is no response after 1.5 seconds of long press, please check if there is any following problem: 1. Is the battery power sufficient? 2. During charging, does the LED light work normally? 3. Is the product too hot?
5 No response when charging	If the instrument does not response to charging, you can remove the trouble as the following steps: 1. Check whether the charging adaptor is abnormal. 2. Is the charging cable and connector damaged? 3. Does the charging LED indicator work?

About Warranty

1. The main body of the product (excluding the free gifts and accessories, such as gel and charging cable) enjoys one-year warranty from the date of purchase.
2. The return and change policy of inFace is in accordance with the service standard of "New Three Guarantees National Policy (China)".
3. If there is any quality problem after opening the box, please contact the after-sales customer service in time.

Limitation of Warranty

1. Failure to use or maintain the product according to the product manual;
2. Machines repaired by non-inFace authorized repair center;
3. Accessories, such as the charging line, free gifts and outer packing are not subject to warranty service;
4. Damage caused by natural disasters and other force majeure;
5. Product failure or damage caused by accident or human factors (including operation error, incorrect insertion and extraction, scratch, disassembly, collision, improper voltage input, etc.).

Notice: Do not immerse the instrument in water at any time. Machine damage caused by human error is not covered by the warranty. If there is any change to the above list, the company's new provisions or the country's updated relevant laws and regulations shall prevail. Please consult the national service hotline 400-007-8388 for the specific warranty terms.



inFace

WARRANTY CARD

Product Number: _____

Product Name: _____

User Name: _____

Phone Number: _____

User Address: _____

Dealer Information: _____

Date of Purchase: _____

Order Number: _____

Trouble Description: _____

(Fill in the warranty card details when applying for the repairing service. The warranty period is one year for the whole product.)

inFace

Product Name: inFace RF Beauty Instrument

Product Model: MS6000

Production Date: _____

Quality Inspection: _____

Shenzhen Migoosmart Co., Ltd.

CERTIFICATION