Operator's Manual should be consulted before using this autoclave.

BASIC OPERATION

1 Power unit ON.

2 Fill with demineralised/distilled water.

NEVER USE TAP WATER

3 Load configurations.

4 Close unit.

5 Run cycle.

A Select a cycle.

B Start a cycle.

C Allow cycle to run.

D End of cycle.

6 Open unit.

7 Unload instruments.

WARM-UP MODE

- STAND-BY LOW – For low to average use: 70°C / 158°F
- STAND-BY HIGH – For high use: 120°C / 248°F
- STAND-BY OFF (default) – For infrequent use. The wait time will be longer (up to 12 minutes from a cold start)

INSTRUCTIONAL VIDEOS

LED COLOUR

WHITE
Unit is idle; ready for operation

RED
Cycle in operation

BLUE
Cycle complete

ORANGE
Cycle stopped or cycle fault condition

LOCK ICON

Drawer is locked; the chamber is still under pressure

Drawer is safe to open

Drawer remains locked; not safe to open

Drawer remains unlocked; unsafe to start cycle

Spare Parts and Accessories	
01-116292S	Air Filter
01-116293S	Drawer Seal, 6L
01-116294S	Tray, 6L
01-116296S	Exhaust Tube
01-116297S	Extended Exhaust Tube
01-116298S	Direct-to-Drain
01-116299S	Water Filling
01-116300S	Waste Bottle
01-116301S	USB Memory Stick
01-116302S	Water Reservoir Cap
SCI-BDSK134V	Bowie-Dick Test Kit
97902001	Helix Test Kit with 100 Chemical Indicators
97902002	Helix Text Kit with 400 Chemical Indicators
01-110282S	Power Cord EU 16A/250V
01-110285S	Power Cord UK IND plug 16A/250V
01-116496	StatMat Silicone Mat STATIM B
01-116497S	Reservoir Chamber Filter
01-116498S	Pouch Rack Kit

IMPORTANT!
Users of this device must be sure to select the most appropriate sterilisation cycle according to the recommendations of the leading infection control authorities and local regulatory guidelines / recommendations.

TOUCHSCREEN

Solid / Unwrapped:
IUSS* cycle for unwrapped solid instruments

Solid / Wrapped:
Wrapped Solid, non-lumened instruments

Universal:
Wrapped and unwrapped solid and hollow instruments, cassettes, textiles

Rubber & Plastic:
Unwrapped and wrapped solid or hollow instruments constructed of metal, rubber and plastic

Helix/BD Test

Vacuum Test

Custom:
For devices that require different cycle parameters than those available. Set according to the instrument manufacturer's instructions.

*Immediate Use Steam Sterilization

MAINTENANCE SCHEDULE

Daily	Wipe the drawer seal .	Monthly or 160 Cycles	Inspect the water reservoir filter (clean or replace).
	Clean external surfaces .		Clean the external fill water tank .
Weekly	Clean the chamber and trays .	Every 6 Months or 1,000 Cycles	Perform all monthly cleaning tasks .
	Drain water reservoir .		Replace bacteriological air filter .
	Disinfect external surfaces .		Replace drawer seal .
	Inspect and clean the chamber filter (or replace if necessary).		
Every 2 Years or 2,000 Cycles	A complete maintenance of the autoclave including testing of the pressure relief valve and the power failure mode (by a SciCan-approved technician) is recommended.		

www.scican.com

95-116554 EU EN R1. STATIM is a trademark of SciCan Ltd. A Coltene Group Company. Manufactured by: SciCan Ltd, 1440 Don Mills Road, Toronto ON M3B 3P9 / CANADA

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