Gebruiksaanwijzing - Elektrische Tandborstel
User Manual - Electric Toothbrush
Manuel d'utilisation - Brosse à Dents Électrique

NL: p.2 - p.7
EN: p.8 - p.13
FR: p.14 - p.20

Model: LG975

CE RoHS

MADE IN CHINA

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▶ Introduction

- Thank you for choosing our LifeGoods product! To ensure proper and safe functioning of the product, please read the user manual carefully before use. Please always follow the instructions for safe usage. If you have any questions about the product or experience any problems, please contact us by email: service@lifegoods.nl.
- The LifeGoods Electric Toothbrush has 5 modes to suit different teeth and gums: White, Clean, Sensitive, Polish, and Massage. Suitable for children, adolescents and elderly people. The Electric Toothbrush is compatible with any power source. The Smart Timer indicates when the 2 minutes of brushing that dentists recommend are over, it also pauses briefly every 30 seconds to indicate how long you have been brushing.

▶ Warning and Safety Instructions

- ⚠ Please ensure that the packaging is kept out of the reach of children! Risk of suffocation!
- ⚠ Make sure to keep this user manual. Read the manual carefully before using your LifeGoods Electric Toothbrush!
- If you have not used the LifeGoods Electric Toothbrush before, it is important that you read and follow the instructions as stated in the user manual. Be aware of the risks and consequences resulting from misuse of the Electric Toothbrush. Misuse of the Electric Toothbrush can result in personal injury, or destruction of the appliance and of the installations.
- Use the Electric Toothbrush only for purposes described in this user manual.
- If the Electric Toothbrush malfunctions, is dropped or (appears to be) damaged in any manner, do NOT use the Electric Toothbrush and contact Customer Service.
- The Electric Toothbrush is designed and intended for private, non-commercial use only.
- Do NOT disconnect the Electric Toothbrush from its power supply by pulling the cable.
- Do NOT carry or pull the device by its cable. The cable should not be kinked, trapped or twisted.
- This device should NOT be used by persons (including children) with physical, sensory or mental disabilities, or persons (including children) with lack of experience with and knowledge of the device, unless they are supervised or have been given instruction concerning the use by a person responsible for their safety.
- Keep the Electric Toothbrush away from children when using and storing.

- People who have a severe periodontal disease or who have had oral surgery should consult a doctor before using the Electric Toothbrush.
- Consult your dentist if excessive bleeding occurs after using this toothbrush or if bleeding continues to occur after 1 week of use. Also consult your dentist if you experience discomfort or pain when you use the Electric Toothbrush.
- The Electric Toothbrush has only been designed for cleaning teeth, gums and tongue. Do NOT use it for any other purpose. Stop using the appliance and contact your doctor if you experience any discomfort or pain.
- Gently place the bristles against the gums and do NOT push too hard.
- Do NOT reconstruct, repair or disassemble the Electric Toothbrush. This can damage the device. For safety reasons, these parts may only be replaced by authorized service personnel. In order to prevent hazards, always send the Electric Toothbrush to the service centre for repairs.
- Make sure your outlet supplies the correct power voltage to avoid damage to the Electric Toothbrush.
- Always check the power cord and plug before using the Electric Toothbrush. Do NOT operate the Electric Toothbrush if the cord or plug are damaged. A damaged cord or plug must be replaced by a qualified person in order to avoid damages and personal injuries.
- Do NOT submerge the Electric Toothbrush in water or any other liquids.
- Keep the Electric Toothbrush away from moisture, heat and in or close to areas where the surrounding temperature is relatively hot (e.g. stoves or other sources of heat).
- Do NOT expose the Electric Toothbrush to direct sunlight.

▶ Receiving Your Shipment


- When receiving your shipment, make sure that the contents are consistent with the packing list. Notify your distributor of any missing items. If the equipment appears to be damaged, file a claim immediately with the carrier and notify your distributor at once, giving a detailed description of any damage. Save the damaged packing container to substantiate your claim.

▶ Contents of Package

- LifeGoods Electric Toothbrush - Model LG975 - EAN 8720195250135

Included Accessories


- 1x User Manual
- 4x Brush Heads
- 1x Protection cap
- 1x EU Adapter

Product Features

Product Overview

1. Brush Head
2. Power Button
3. White Mode
4. Clean Mode
5. Sensitive Mode
6. Polish Mode
7. Massage Mode
8. Charging Light
9. Charging Port

Technical Data

Size	25,60 x 25,50 x 227 mm
Weight	70,40 gr
Colour	Black
Material	ABS / TPE
Cable Length	100 cm
Input	DC5V,1.5W ,0.3A
Output	DC3.7V,1W
Battery Type	Lithium, DC3.7V 260 mAh
Power Consumption	1.0 W

Bristle Type	Dupont Nylon
Number of modes	5 (White, Clean, Sensitive, Polish, and Massage)
IPX value	IPX7
Working Temperature	0-50°C
Charging Time	± 4 hours (lasts 3 weeks)
Strokes per minute	40,000
Automatic shut-off	Yes, after 2 minutes

Other Features

- The bristles change colour to indicate when you have to change the bristle;
- Hygienic protection cap included;
- With a helpful Smart Timer.

Operating Instructions

 The Electric Toothbrush has only been designed for cleaning teeth, gums and tongue. Do NOT use it for any other purpose. Stop using the appliance and contact your doctor if you experience any discomfort or pain.

 Gently place the bristles against the gums and do NOT push too hard.

How to Use



1. Attach the brush head and make sure it is locked in place.
2. Press the Power button to choose your preferred mode. You can shift between the 5 modes “White, Clean, Sensitive, Polish, Massage” by pressing the Power button continuously.
3. Wet the bristles and apply the required amount of toothpaste.
4. Press the Power button to turn the Electric Toothbrush on.
5. Gently move the brush head across the front of your teeth, the back of your back teeth and your molars. Gently keep the bristles placed on the teeth and in the gumline. Brush your teeth with a small back and forth motion so the bristles reach between the teeth. Do NOT press the bristles too hard on the same area for an extended period of time.
6. The recommended amount of time for brushing is 2 minutes. After every 30 seconds, there is an interval pause to remind you how long you have been brushing.
7. To make sure you brush evenly throughout the mouth, divide the mouth into 4 sections. Move to the next section when you feel the interval pause. The Electric Toothbrush will automatically power off after 2 minutes.

8. To avoid splashing toothpaste, turn the Electric Toothbrush on after you put it against your teeth, and only remove it from your mouth when the toothbrush is turned off.

Changing the Brush Head

- **Install:** Slide the brush head down to the handle. Turn clockwise until A meets T.
- **Detach:** Turn brush head counter clockwise until A meets B. Gently pull the brush head upwards.

How to Charge

1. When the battery power is low, the charging indicator light will flash 10 times to remind you to charge your Electric Toothbrush.
2. Use the USB cable to connect the Electric Toothbrush to a power source to start charging.
3. When the battery is fully charged, the charging indicator light will turn off and the charging unit will shut off automatically to avoid overcharging.
4. When you use the Electric Toothbrush for the first time or have not used the Electric Toothbrush for an extended period of time, the charging indicator light will not light up for the first 5 minutes.

Cleaning and Maintenance

- Remove the brush head and rinse the bristles and the metal shaft area with warm water. Clean the bristles after every use, and clean the metal shaft area at least once a week.
- Wipe the entire surface of the handle with a damp cloth.

Note: Do NOT use rubbing alcohol, vinegar or bleach to clean the handle as this may cause discoloration or might otherwise damage the Electrical Toothbrush.


Storage

- If you are not going to use your LifeGoods Electric Toothbrush for an extended period of time, unplug it from its powersource, clean it and store it in a cool and dry place away from direct sunlight.


Disposal - Recycle

- Contribute to a cleaner living environment! Do not dispose of your LifeGoods Electric Toothbrush with the household waste. For proper recycling, please take these products to a designated collection point. Please contact your local authority or your household waste disposal service for further details of your nearest designated collection point.

Service and Warranty

- If you need service or information regarding your product, please contact the LifeGoods Customer Support at service@lifegoods.nl
- LifeGoods provides a 2 year warranty on its products. To obtain service during the warranty period, the product needs to be returned with a proof of purchase. Product defects have to be reported within 2 years from the date of purchase. This warranty does not cover: damages caused by misuse or incorrect repair; parts subject to wear; defects that the customer was aware of at the time of purchase; damages or defects caused by customer neglect; damages or defects caused by third parties.

Disclaimer

- Changes reserved; specifications are subject to change without stating grounds for doing so.

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