ARDESTO

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IMPORTANT SAFEGUARDS

Please read the safety instructions before using this product. The safety instructions provided here are to help you use this product safely and correctly, and also to prevent you from injuring yourself or others.

READ ALL INSTRUCTIONS BEFORE USING

WARNINGS

To reduce the risk of burns, fire, or injury to persons:

- This device is not a toy. Children should not play with the device. Cleaning and customer service should not be performed by children without supervision.
- Use this product only for its intended use as described in this manual. This product is designed to clean your teeth, gums only. Do not use any attachments or brush heads other than those recommended by the manufacturer.
- The toothbrush is a personal care appliance and is not intended for use on multiple patients in a dental office or institution.
- If your toothbrush includes peroxide, baking soda or other bicarbonate (common in whitening toothpastes), thoroughly clean your brush head and handle to avoid cracking the plastic.
- Keep the product away from heated surfaces.
- If you experience discomfort or pain in your mouth, discontinue use of this product and contact a dentist.
- Do not clean the product in the dishwasher.
- Do not use outdoors or operate where aerosol (spray) products are being used or where oxygen is being administered.
- Never drop or insert any object into any opening.
- For household use, indoor use only.
- If long time don't use, please remove the batteries, so as not to cause danger or damage to the products such as batteries leakage.

MEDICAL-WARNINGS

- If you have had oral surgery or gum surgery within the past two months, consult your dental provider before using the product.
- Consult your dentist if excessive bleeding occurs while using this product or minor bleeding continues to occur after one week of use.
- Consult your physician or implant device manufacturer prior to using this product if you have questions concerning the combined use of this product with a pacemaker or other implanted medical device(s).
- Consult your physician prior to using the toothbrush if you have any medical concerns.

Production date and device serial number are indicated on the packaging.

SAVE THESE INSTRUCTIONS

EQUIPMENT

Electric toothbrush - 1 pc.


Brush head - 3 pcs.

Protective cap - 1 pc.

Operating instructions (contains a warranty card) - 1 pc.

CONFIGURATION

The Electric Toothbrush For Kids is designed for children more than 5 years old. It is not intended for adult use.

For children under age 8, dental professionals recommend that tooth brushing is supervised by a parent or adult caregiver.

GETTING STARTED

Install the batteries. (Fig. 1).

WARNING! Batteries should only be installed by adults. Batteries cannot be charged. Rechargeable batteries must be charged before use. Rechargeable batteries should only be charged under the supervision of an adult. For best performance, use the recommended batteries. Observe the polarity when installing the batteries. Power terminals must not be short-circuited.


Fig. 1

HOW TO USE

- ① Slide the brush head on the toothbrush according to the mask, insert it into the drive shaft catch and turn it clockwise, until it locks into place on handle. (Fig. 2).

- The prevents the brush head from becoming detached due to vibration.

- ② Put toothpaste on the brush head.

- ③ Place the brush head in your mouth.

- ④ Lightly press the on button to switch on the toothbrush and start to brush your teeth. (Fig. 3). The toothbrush manual switches off at end of the brushing cycle.

- Brush evenly throughout the mouth in four sections. (Fig. 4).

1. Outside upper; 2. Outside lower;
3. Inside upper; 4. Inside lower;

- ⑤ According to brushing time and sections, guide the brush head slowly from teeth to teeth.

- Hold the brush head in place for a few seconds before moving on to the next tooth.

- Do not press too hard or scrub, simply let the brush do all the work.

- ⑥ During brushing, lightly press off button to terminate a brush. (Fig. 5).


Fig. 2


Fig. 3


Fig. 4

NOTE:

- A brushing time of 2 minutes is generally long enough to clean your teeth properly. If you want you can brush longer than 2 minutes for an even more thorough clean.
- When you first use you electric toothbrush you may experience slight bleeding. This should stop after a few days when your gums get used to the toothbrush. If the bleeding persists, Please consult your dentist.


Fig. 5

MAINTENANCE

- ① Turn the brush head counter-clockwise, remove the toothbrush head from toothbrush handle. (Fig. 6).
- ② Turn on the toothbrush and rinse the drive shaft under the tap to remove toothpaste residues. (Fig. 7).
- ③ Turn off the toothbrush, detach the brush head and rinse the inside of the brush head (max. 140F). (Fig. 8).
- ④ Turn on the toothbrush and rinse the brush head under a tap. (Fig. 9).
- ⑤ Toothbrush head put back on the handle. (Fig. 2).
- ⑥ Clean the toothbrush with a soft moist cloth or paper towel. (Fig. 10).

When the toothbrush is not used for a long time, you should be take out the batteries, and place the product into the packaging box.


Fig. 6


Fig. 7


Fig. 8

NOTE:

- Replace the brush head every 3 months for the best brushing effect.
- Also replace the brush head when the bristles start to bend outwards.
- Any other servicing should be performed by an authorized service representative or that the product has no user serviceable parts.
- Do not mix old and new batteries.
- Do not mix different types of batteries.


Fig. 9


Fig. 10

TROUBLESHOOTING

Problem	Possible Cause	Solution
The toothbrush doesn't work	No batteries inserted	Insert batteries in the unit
	The batteries in wrong direction	Re-insert batteries in correct direction
	The batteries are exhausted	Replace with new batteries
Weak or diminished oscillations	The toothbrush is pressed hard against your teeth	Place the brush head gently against your teeth

CORRECT DISPOSAL OF THIS PRODUCT

This marking indicates that this product should not be disposed with other household wastes. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

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Warranty conditions:

The plant guarantees normal operation of the product within 12 months from the date of the sale, provided that the consumer complies with the rules of operation and care provided by this Manual.

Service life of the product is 3 years. In order to avoid any misunderstanding we kindly ask you to carefully study the Owner's Manual, the terms of warranty, and check the correctness of filling the Warranty Card. The Warranty Card is valid only if the following information is correctly and clearly stated: model, serial number of the product, date of sale, clear seals of the company-seller, signature of the buyer. The model and serial number of the product must correspond to those specified in the Warranty Card. If these terms are violated, and if the information specified in the Warranty Card is changed, erased or rewritten, the Warranty Card is considered invalid. Setup and installation (assembly, connection, etc.) of the appliance are described in the attached documentation; they can be carried out both by the user, and by specialists having the necessary qualification sent by USC or selling companies (on a paid basis). In this case, the person (the institution) installing the appliance shall be responsible for the correctness and quality of installation (setup). Please pay attention to the importance of the correct installation of the appliance, both for its reliable operation, and for obtaining warranty and free service. Require the installation specialist to enter all the necessary information about the installation in the Warranty Card. In case of failure of the appliance during the warranty period by the fault of the manufacturer, the owner is entitled to free warranty repair upon presentation of the correctly filled Warranty Card together with the appliance in the factory configuration packed in the original packaging to the warranty workshop or the place of its purchase. Satisfaction of the customer's claim due to the fault of the manufacturer shall be made in accordance with the law «On Protection of Consumers' Rights». In the case of warranty repair the warranty period shall be extended for the period of repair and shipment.

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WARRANTY CARD

ARDESTO

Dear Buyer! Congratulations on your purchase of the ARDESTO brand appliance, which was designed and manufactured in accordance with the highest quality standards, and we thank you for choosing this particular instrument.

We ask you to keep the coupon during the warranty period. When purchasing a product, require a full warranty card.

1. Warranty service is carried out only if there is a correctly and clearly filled original warranty card, which indicates: product model, date of sale, serial number, warranty service period, and the seller's seal. *
2. The service life of household appliances is 3 years.
3. The product is intended for consumer use. When using the product in commercial activities, the seller / manufacturer does not bear warranty obligations, after-sales service is performed on a paid basis.
4. Warranty repair is carried out within the period specified in the warranty card for the product in an authorized service center on the conditions and terms determined by applicable law.
5. The product is withdrawn from the warranty in case of violation by the consumer of the operating rules set forth in the instruction manual.
6. The product is removed from warranty service in the following cases:
 - misuse and non-consumer use;
 - mechanical damage;
 - damage caused by the ingress of foreign objects, substances, liquids, insects;
 - damage caused by natural disasters (rain, wind, lightning, etc.), fire, domestic factors (excessive humidity, dust, aggressive environment, etc.)
 - damage caused by non-compliance of power and cable network parameters with state standards and other similar factors;
 - when operating equipment in the power supply network with a missing single ground loop;
 - in case of violation of seals installed on the product;
 - lack of serial number of the device, or inability to identify it.
7. The warranty does not cover consumables and accessories, as well as filters, shelves, drawers, stands and other containers for storing products.
8. The warranty period is 12 months from the date of sale.

* Tear-off maintenance tickets are provided by an authorized service center.

The completeness of the product is checked. I have read the terms of the warranty service, no complaints.

Customer Signature _____

Elit service :

**3, Petre Kavtaradze str., Tbilisi, 94, Gorgiladze str., Batumi
71, Davit Guramishvili str., Gori, 73 Kostava str., Akhaltsikhe**

Product Information

Product

Model

Serial number

Seller Information

Name of trade organization

Address

Date of sale

Seller stamp

Warranty coupons

Ticket number 3

Seller stamp

date of the application

Cause
of damage

Date
of execution

Ticket number 2

Seller stamp

date of the application

Cause
of damage

Date
of execution

Ticket number 1

Seller stamp

date of the application

Cause
of damage

Date
of execution