·IDUTION·

ORAL IRRIGATOR

DT-FC2660

USER MANUAL

ENGLISH	4
DEUTSCH	9
ČESKY	14
SLOVENSKY	19
MAGYAR	24

SAFETY INFORMATION

Please read this manual carefully before use and keep it for future reference. The manufacturer declines any liability for damage caused by improper use or failure to follow the instructions in this manual. Improper use or non-compliance with the instructions will void the warranty.


- Oral irrigators can be used by children aged 8 and over and people with reduced physical or mental abilities or lack of experience and knowledge if they are under supervision.
- Cleaning and user maintenance must not be performed by children without supervision.
- Children must not play with the device.
- Oral irrigators must only be used with cold or hot water or with special solutions as defined in the manufacturer's instructions.
- If the power cord, plug, or any other part of the product malfunctions, or if any part of the product has been dropped or otherwise damaged, repair may only be performed by a qualified electrician. Improper repairs can put users at risk of injury.
- The device must be positioned so that there is no risk of it falling into water or on the ground. Do not use the oral irrigator in the bath or shower.
- Slight gum bleeding may occur the first time you use the oral irrigator. However, if the bleeding is heavy or does not stop within a few days, stop using the device and consult a dentist.
- If you have recently had dental or gum surgery, suffered from long-term dental problems or have undergone extensive procedures in the last 2 months, or if you suspect dental problems (fillings, crowns, etc.), consult your dentist before using the device.
- The device meets the safety standards for electromagnetic devices. If you have a pacemaker or other implanted device, consult your physician or the manufacturer of the implanted device before use.
- Do not aim the jet under your tongue, ear, nose, eyes or other sensitive areas, as this may cause injury.
- This device is intended for personal care, not for use by multiple patients in a dental office or institution.

TECHNICAL SPECIFICATIONS

Dimensions	71.5 × 63 × 220 mm
Weight	235 g
Input	5 V / 1 A
Battery	Li-Ion, 3,7 V 1400 mAh
Protection level	IPX7
Water tank volume	260 ml
Maximum flow	220 ml/min
Battery life	40 minutes

OVERVIEW

- | | |
|--|-----------------------|
| 1. Nozzle | 5. Mode change button |
| 2. Nozzle ejection button | 6. Water tank |
| 3. Main switch | 7. Charging connector |
| 4. Display showing irrigator mode and battery capacity | 8. Filling opening |

Package components not shown: power cable, replacement nozzle

USING ORAL IRRIGATOR

Filling with water and starting

- Before using for the first time, allow the oral irrigator to be fully charged using the included cable. Keep in mind that the oral irrigator cannot be used while charging.
- Open the cover on the back of the water tank and fill the tank with an adequate amount of water.
- To start the oral irrigator, press the main switch, press the switch again to switch off the irrigator.

Switching oral irrigator modes

- After starting by pressing the switch, the oral irrigator will start in the mode in which it was switched off during the previous use.
- To switch the oral irrigator mode, press the mode change button.
- The current mode will be shown on the display by a pictogram and a name.
- Oral irrigator offers 4 modes: **STRONG**, **MEDIUM**, **SOFT** and **PULSE**. **STRONG** mode has the strongest intensity, **MEDIUM** medium and **SOFT** the gentlest. In the **PULSE** mode, the intensity of the jet will change between high and low performance.

Cleaning tips

- Adjust the position of the irrigator nozzle so that it points at the gum line at an angle of approx. 90°.
- Keep your mouth closed while cleaning the interdental spaces to prevent water splashing around.
- Clean the interdental spaces in the order shown in the picture. Start at the back and continue along the gums forward, to the back teeth on the opposite side of the oral cavity.

CLEANING AND MAINTENANCE

- Always pour all residual water out from the tank after each use of the oral irrigator.
- If you will not be using the oral irrigator for a long time, clean all parts and store them in a cool, dry place out of direct sunlight.

Charging

The oral irrigator is equipped with a rechargeable battery, which can be easily recharged using the included cable. The remaining battery capacity is shown on the display after switching on.

- Plug the USB cable into a suitable charger with a USB connector.
- Connect the other end of the cable with the connector to the connector of the oral irrigator.
- Charging takes about 4 hours.

TROUBLESHOOTING GUIDE

Problem	Possible cause and solution
Water leaks between the irrigator attachment and the body of the appliance	The attachment is not inserted correctly. You can hear the correct insertion of the attachment by an audible snap.
The water pressure is insufficient.	The battery is probably empty. Charge the battery with the included cable.
Water leaks from the tank	The seal is probably damaged. Contact a service centre or vendor.
The oral irrigator cannot be switched on	The battery is probably empty. Charge the battery with the included cable.
No water flows out of the nozzle	The irrigator nozzle seal is probably damaged. Contact a service centre or vendor.

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ENGLISH

The use of the WEEE Symbol indicates that this product should not be treated as household waste. By ensuring this product is disposed of correctly, you will help protect the environment. For more detailed information about the recycling of this product, please contact your local authority, your household waste disposal service provider, or the shop where you purchased the item.

This product meets all of the related basic EU regulation requirements. The EU declaration of conformity is available on www.alzashop.com/doc.

This product sold in the European Union meet the requirements of Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS).

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