SVENSKA

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PROCARE®

METAL FINGER SPLINTS/
FINGER PROTECTORS/FINGER COTS/
FROG SPLINTS/LEWIN STERN THUMB/
FINGER SPLINT

SPLINT DI METALLO PER DITA/PROTEZIONI PER DITA/COPRIDITA IN PELLE/SPLINT GONFIABILI/SPLINT LEWIN STERN PER POLLICE/DITA
FÉRULAS METÁLICAS PARA LOS DEDOS/PROTECTORES PARA LOS DEDOS/DEDILES/FÉRULAS EN FORMA DE RANA/FÉRULA PARA PULGAR/DEDOS LEWIN STERN

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
ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE:

The intended user should be a licensed medical professional, the patient or the patient's caregiver. The user should be able to read, understand and be physically capable of performing the directions, warnings and cautions in the information for use.

INDICATIONS:

For immobilization of the fingers/thumb following soft tissue injuries, fractures or dislocations. Providing immobilization or controlled movement of the limb or body segment.

APPLICATION INSTRUCTIONS:

A) Place splints on finger/thumb in desired position and secure with tape.

Lewin Stern Splint

B) Secure wrist strap.

USE AND CARE:

The splints are malleable and may be contoured for individual requirements.

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

WARNINGS AND PRECAUTIONS:

- If you have skin disorders/injuries in the relevant part of the body, the product should be fitted and worn only after consultation with a medical professional.
 - Do not use over open wounds.
 - If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact your doctor immediately.
 - Do not use this device if it was damaged and/or packaging has been opened.
- NOTE:** Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

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