



Quick Start Guide

personal microderm elite pro

CLINICAL-GRADE EXFOLIATION WITH VACUUM SUCTION

CAUTION: Please read this manual carefully before using

GETTING STARTED

Customize Your Treatment

The Personal Microderm Elite Pro is ideal for all skin types.

1. TRAVEL LOCK

Device may come with travel lock activated. To disable, hold the power and + buttons simultaneously for 3 seconds. The lights on the front will flash once for unlocking and twice for locking.

2. SELECT CAP



Small Luxe Cap
for face or smaller areas



Large Luxe Cap
for body



Blackhead Cap
for face

3. SELECT DISC

**Do not use disc with blackhead cap.*




Small Disc
for face or smaller areas



Large Disc
for body

From ultra sensitive
to very intense.



*Note: not all discs are included in the set.
Learn more at pmdbeauty.com/discs*

4. SELECT MODE

- Low Speed & Low Suction: **Sensitive Treatment**
- Moderate Speed & Moderate Suction: **Moderate Treatment**
- High Speed & High Suction: **Intense Treatment**

Replace Discs Every 3-4 Uses

After consistent treatments, skin may become more tolerant. Gradually increasing disc intensity is recommended.



HOW-TO USE

Watch our how-to video



TIPS & TRICKS

- 1 Start with white training disc
- 2 Glide quickly along skin
- 3 Do not apply pressure
- 4 Do not pause in one area
- 5 Only one pass per area
- 6 Press down disc firmly on post




CHARGING



Charging Time:
3 hours

Device will flash 3 times
for low battery signal.



Plug the USB end of the
charging cord into any
outlet adapter. Insert
the opposite end of the
charging cord into the
base. Place device on
charger.



When device is connected
to the charging base it will
vibrate twice. Flashing
battery light indicates
charging status. When
fully charged, lights stop
flashing.

FAQ's

When can I expect to see results?

We promise that after your first treatment you will have softer, smoother skin. After eight to twelve weeks of consistent use, individuals may see a reduced appearance of fine lines and wrinkles, dark spots, pore size, and uneven skin tone and texture.

How often should I be doing a treatment?

The Personal Microderm Elite Pro should only be used once a week. Microdermabrasion removes the top layer of dead, dull skin, and your body needs one week to rejuvenate. Be consistent, you will see results.

CLEANING & MAINTENANCE

Rinse the cap and filter between each use with warm water and soap. Be sure to replace the disc every 3-4 uses and replace the filter every 3-4 months. Wipe down the body of the device as needed and store your Personal Microderm Elite Pro in a dry place. Expect to charge the device every 4-6 months with regular use.

Caution and Warnings:

The effects of the Personal Microderm Elite Pro are not always immediately visible. Not following instructions including: multiple passes or pausing in one spot may lead to superficial red marks or even light scabbing. If red marks or light scabbing occur, do not be alarmed—they will disappear. Review the instructions again and wait 7 days until the next treatment. Do not use any active or peeling skincare products 48 hours before or after using the Personal Microderm Elite Pro—your skin may be too sensitive. If you have any questions before beginning a treatment, please reach out to our experts at 888-445-4532.

#BRILLIANTCONFIDENCE

To learn more visit: pmdbeauty.com/elite

FOLLOW US

YouTube
@pmd

Instagram
@pmdbeauty

Pinterest
@pmdbeauty

