


ANJOU



MODEL: AJ-PCA028

PORTABLE ORAL IRRIGATOR User Guide

Thank you for choosing Anjou product. Please register your product at our website www.anjou.com for further support.

IMPORTANT

Safeguards

All information in this section is provided to instruct you on the correct and safe operation of this device so that you can avoid injuries to yourself and others as well as damage to property.

Dangers


- Toddlers, infants, people who are unable to operate should not be allowed to use this appliance. Please keep it out of reach of children.
- Never alter the appliance. Also, do not try to disassemble or repair it.
- Do not place or store the product where it can fall or be pulled into tub or sink.

Warnings


- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
- The appliance is only to be used with the power supply unit provided with the appliance.
- The battery must be removed from the appliance before it is scrapped.
- The appliance must be disconnected from the supply mains when removing the battery.
- The battery is to be disposed of safely.

- Regarding how to remove the battery, please see the "Battery Maintenance and Disposal" section.
- Do not fill the tank with water warmer than 40 °C / 104 °F.
- If the cable is damaged, it shall be replaced by the manufacturer or its service agent or a similarly qualified person in order to avoid a hazard.
- Keep the cord away from heated surfaces.
- Please do not use it if the nozzle is lost.
- Do not use for any purpose other than oral cleaning.
- Keep nozzle away from children. The nozzle is small enough to be swallowed and can be lodged in the throat.
- People who experience tooth or gum pain may not be able to use this device and should consult their dentist.
- Do not press the nozzle too firmly against your teeth or gums. Doing so may injure your gums.

PACKAGE CONTENTS



PRODUCT DIAGRAM



- | | |
|-------------------------|----------------------------|
| ① Nozzle | ⑧ Battery indicator |
| ② Nozzle Release Button | ⑨ Charging Port |
| ③ Power Button | ⑩ Reservoir |
| ④ Mode Button | ⑪ Air Inlet |
| ⑤ Normal Mode Indicator | ⑫ Water Inlet |
| ⑥ Soft Mode Indicator | ⑬ Detachable Reservoir Lid |
| ⑦ Pulse Mode Indicator | |


SPECIFICATIONS

Model	AJ-PCA028
Input	DC 5V --- 1A
Power	5W
Charging Time	Approx. 4 hours
Battery Life	2 weeks of lifetime between charges
Reservoir Capacity	300 ml / 10.0 fl. oz
Built-in Battery Capacity	1400mAh
Working Cycle	2 minutes
Dimensions	8.5 x 7.2 x 22.6 x 24.9cm / 3.3 x 2.8 x 9.8in
Weight	249g / 8.8oz

GET STARTED

Charging

1. Connect the charging port with the included charging cable.
2. Connect it to any USB output port, and charging will start automatically.
3. The charging indicator will flash blue during charging and flash green once fully charged.




Note:

- Please fully charge the unit before first use, which will take around 4 hours.
- It is not recommended to use the unit when charging.

HOW TO USE

Mounting and Removing the Nozzle




1. Insert the nozzle into the hole at the top of product. Push the nozzle in vertically, make sure it is well aligned with the body to avoid water leakage.
2. To remove the nozzle, press the Nozzle Release Button and pull away from the product.

Note:

- When the nozzle is mounted, rotate it as per your preference.
- Change the nozzle every six months or when the nozzle can't fit tightly into the body.

Adding Water



Using one of the following ways to add water:

1. To add water, hold the product horizontally, open the reservoir lid, and add water through the inlet, lock the lid into place.
2. Cap the lid, remove the reservoir, and add water from the top. Install the reservoir back on.

Note:

- If the lid falls off, reinstall it as instructed (refer to diagram 3).
- You can add warm water into the reservoir in winter (refer to diagram 4).



www.anjou.com

NORTH AMERICA

E-mail: support@anjou.com






MANUFACTURER

Manufacturer: Shenzhen NearbyExpress Technology

Development Co., Ltd.

Address: Floor 7, Building E, Galaxy World Phase II, Shenzhen, China

Learn more about the EU Declaration of Conformity:

<http://customer.anjou.com/media/downloads/AJ-PCA028-CE-Cert.pdf>**Power On / Off**

1. Press the Mode Button to cycle through Normal, Soft, and Pulse Mode and the corresponding indicators will light up.
2. Lean over the sink, place the nozzle in your mouth, aim the nozzle towards your teeth, and then turn the unit on.





Note:

- It is recommended to choose the Soft mode for the first use, then adjust the modes as per your preference after 2-5 uses.
- The unit will remember the mode setting of last use when turned on again. Please adjust it as previously instructed if needed.

3. Direct the nozzle at a 90-degree angle to your gum line with your mouth partially closed to avoid splashing, but allow water to flow from your mouth into the sink.
4. For the best results, start with the molar area (the back teeth) and move along the gum line toward the front teeth, pausing briefly between teeth.
5. Continue until all areas around and between the teeth have been cleaned.

Note:

- The irrigator stops automatically after the 2 minute working cycle is done(refer to diagram 6).
- Please keep the nozzle inside your mouth if you want to turn off the unit during a working cycle to avoid splash.

CLEANING

1. Remove the nozzle and rinse it under tap water to clear all the debris.
2. Remove the reservoir and rinse it under tap water.
3. Brush the filter gently with soft brush.
4. Wipe the body with wet cloth.

Note:

- Make sure to unplug the unit before cleaning.
- Do not wash the unit in dishwasher.

BATTERY MAINTENANCE DISPOSAL**Battery Maintenance**

If the Oral Irrigator is to be stored for an extended period of time, charge it for 4 hours before storing.

Battery Disposal

The unit contains a Nickel Metal Hydride battery. At the end of the battery life, please recycle the unit at your nearest recycle center per local or state requirements. The battery is not replaceable.

CAUTION

- Check the charging cord for damage before the first use and during the service life of the product.
- Do not direct the nozzle under the tongue, into the ear, nose, or other delicate areas, as this product is capable of producing pressure that may cause serious damage.
- Do not plug this product into a power source that is different from the specifications specified on the device or charger.
- Use this product only as indicated in these instructions as recommended by a dental professional.
- This appliance should not be used by children or individuals with certain disabilities unless under supervision or instruction by another person responsible for their safety.
- Bleeding may happen during the initial use due to the effective cleaning of your gingival sulcus.

TROUBLESHOOTING

Problem	Cause	Solution
Water pressure too strong	Not used to the water pressure for the initial use	Select the Soft mode for the initial use and adjust to other modes after 2-3 uses
Water pressure too weak	Low voltage	Please fully charge the unit before use
Gum bleeding	1. Effective cleaning of the gingival sulcus during the initial use 2. Improper cleaning methods	Clean your teeth as instructed in this manual
No water or water leakage from the nozzle	Did not install the nozzle properly	Please reinstall the nozzle and make sure it is fit tightly into the body
Reservoir lid falls off	The reservoir lid is detachable	Install it as instructed in the "Adding Water" Section
Fail to charge the unit	1. Debris on the charging port 2. Battery malfunction	1. Clean the charging port regularly with Q-tip 2. If charging still doesn't start, consult Anjou after sales support team

**WEEE Compliance**

Correct Disposal of This Product (Waste Electrical & Electronic Equipment)(Applicable in countries with separate collection systems)
This marking on the product, accessories or literature indicates that the product and its electronic accessories should not be disposed of with other household waste at the end of their working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate these items from other types of waste and recycle them responsibly to promote the sustainable reuse of material resources.
Household users should contact either the retailer where they purchased this product, or their local government office, for details of where and how they can take these items for environmentally safe recycling. Business users should contact their supplier and check the terms and conditions of the purchase contract. This product and its electronic accessories should not be mixed with other commercial wastes for disposal.

FCC Compliance

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:
(1) this device may not cause harmful interference, and (2) this device must withstand any interference received, including interference that may cause undesired operation.