



5 Fr. **2** Lumen **40** cm catheter length **.018** inch dia. spring-wire guide



Pressure Injectable Two-Lumen PICC

Contents:

- | | | | |
|----|---|-----|--|
| 1: | Two Lumen Catheter: 5 Fr. (1.81 mm OD) x 40 cm, Pressure Injectable, T-Port Connector, Blue FlexTip [®] , Contamination Guard and Placement Wire | 1: | Safety Scalpel: #11 |
| 1: | GlideThru [™] Peel-Away Sheath: 5 Fr. x 2-3/4" (7 cm) Radiopaque over 5 Fr. Dilator | 1: | Patient ID Card |
| 1: | Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer | 1: | Sterile Procedure Sign |
| 1: | Injection Needle: Eclipse [™] 22 Ga. x 1-5/8" (4 cm) | 1: | Patient Information Booklet |
| 1: | Injection Needle: Eclipse [™] 25 Ga. x 1" (2.54 cm) RW | 1: | Chart Sticker |
| 1: | Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW | 1: | PICC Line Insertion Checklist |
| 1: | Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip | 1: | Paper Tape Measure |
| 2: | Syringe: 10 mL Luer-Lock | 1: | Tourniquet |
| 1: | Syringe: 3 mL Luer-Lock | 5: | Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 1: | SecondSite [™] Adjustable Hub: Fastener | 1: | Dressing: Tegaderm [®] I.V. 4" x 6" (10 cm x 15.5 cm) |
| 1: | SecondSite [™] Adjustable Hub: Catheter Clamp | 1: | Dressing: GRIP-LOK [®] Catheter Fixation Device |
| 1: | SharpsAway [®] Disposal Cup | 10: | Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1: | SharpsAway [®] II Locking Disposal Cup | 1: | Surgical Apparel: Cap |
| 1: | Catheter Trimmer | 1: | Surgical Apparel: Impervious Gown |
| 1: | Towel: 24" x 36" (61 cm x 91 cm) | 1: | Dressing: Steri-Strip [®] 1/2" x 4" (1.2 cm x 10 cm) |
| 1: | Drape: 60" x 76" (152 cm x 193 cm) | 1: | Surgical Apparel: Mask |
| 1: | Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration | 1: | Surgical Apparel: Mask with Shield |
| 1: | Towel | 1: | Transducer Cover: 4-12/16" x 95" (12 cm x 244 cm) |
| | | 2: | Valve: MicroCLAVE [®] Neutral Displacement Connector |
| | | 1: | HemoHopper [®] Fluid Receptacle |
- ¹A trademark of Becton, Dickinson and Company.

²A registered trademark of 3M Company.
³A registered trademark of TIDI.
⁴A registered trademark of ICU Medical, Inc.
All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU


Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

Contraindications: The Pressure Injectable PICC is contraindicated wherever there is presence of device related infections or presence of thrombosis in the intended insertion vessel or catheter pathway. Clinical assessment of the patient must be completed to ensure no contraindications exist.

 Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	Pump Flow Rate†† (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (18 Ga.)	0.36	416	3500	4
Proximal (18 Ga.)	0.38	475	4090	4

* Priming volumes are approximate and are done without accessories.
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.
‡ Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.
** Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.



REF EU-24052-HPMSB

5 Fr. **2** Lumen **40** cm catheter length **.018** inch dia. spring-wire guide

