

SKG


T5-1



Scan the QR code to get a quick start

User's Guide

1. Pasting the pain-relief patch




Open the packaging, tear off the white protective film, paste it onto the massage patch and then press it for firmly pasting



Tear off the transparent protective film, then it is ready to use

2. Pain Relief Function of Massager



Single frequency mode of high frequency wave, no stimulation, deep relief of cervical pain and alleviating dowager's hump




The combination of middle frequency/high frequency wave and hot compress has more obvious effect on relieving neck pain

3. Optimal Massage Areas and Sore Spots


The blue area is the best massage area

The optimal sore spots are shown in the picture:

- ① Sore spot of splenius capitis and cervicis
- ② Sore spots of trapezius
- ③ Sore spots of musculus levator scapulae



Just turn it on for a good experience!



Instructions of Power Supply and Charging

Battery level ≥25%	The blue indicator light is normally on, indicating normal battery level
Battery level ≤25%	The orange indicator light is normally on, indicating low battery level
Charging status	The orange indicator light flashes
Fully charged	The blue indicator light is normally on

Technical Parameters

Product Name	Shoulder and Neck Massager	Model	T5-1
Battery Capacity	500mAh	Rated Charging Voltage	5V
Executive Standard	GB4706.1-2005 GB4706.10-2008	Rated Input Power	5W

* The model approval code of the product is displayed in the "product nameplate by laser carving".

Innovative T5-1 Deeply Relieves Muscular Pain



Troubleshooting

Fault	Possible Cause/ Solution
Medium frequency mode enabling failure after startup	1. The two massage patches are not pasted to the skin at the same time. 2. Check whether the switch of the instrument is pressed.
Electric shock feeling up on contacting the massage patch in medium frequency mode	It is a normal phenomenon. The working principle of the massager is that two massage patches form an electrical circuit through the human body. Please make sure that the medium frequency mode is not turned on before moving.
Charging failure or charging indicator light not on	Please check whether the Type-C cable is properly connected to the instrument and the power adapter of the mobile phone.

Notes for Product Scrap

Please recycle the rechargeable battery (lithium ion battery)!

	Danger The battery of this product is special rechargeable battery and shall not be used in other products; besides, it cannot be charged after it is taken out the instrument. - Do not throw it into the fire, resulting in overheating - Do not puncture it with nails, get it to suffer from gravity impact, or disassemble it - Do not connect the positive and negative terminals of the battery directly to metal - Do not store it together with necklaces or hair clips - Do not charge, use or put it in the places in high temperature such as fireside or in the blazing sun - The situations above may cause heating, fire and breakage.
	Warning Please keep the rechargeable battery out of reach of infants to prevent it from being swallowed by mistake. - If the battery is accidentally swallowed, seek medical advice immediately. When the battery is leaking, do not touch the battery with bare hands but deal with it with the following methods. - Do not rub your eyes with your hands when battery electrolyte enters your eyes in case of blindness. Rinse your eyes with clean water immediately and seek medical advice. - When your body or clothing is contaminated by the battery electrolyte, it may cause skin irritation or injury. Please rinse with water. If you still feel discomfort, seek medical advice immediately. - The product can only be disassembled when it is discarded. Failure to do so may result in fire, electric shock, injury or other accidents. If you need to repair the product (except for battery replacement), please consult the shop where you purchased the product.

The lithium ion battery of this product is a valuable resource that can be recycled. Please take out the lithium ion battery and hand it to a professional institution that can recycle and regenerate the rechargeable battery when the product is discarded after use.

- The battery must be removed from the instrument before the instrument is discarded;
- The instrument must be powered off before removing the battery;
- The battery should be disposed of safely.

Removal steps of battery:

- Open the shell with a tool;
- Cut off the battery lead with scissors;
- Take out the battery;

Precautions

Warning

- The product must work with the pain relief patch, otherwise you may get electric shock.
- High frequency waves from the massager may react when the device is close to the store's safety gate.
- If the patient uses the high-frequency surgical equipment and the massager on the body at the same time, burns may occur in the position of the massager patch and the massager may be damaged.
- The output of the massager may be unstable if the product is used near the shortwave or microwave therapy equipment (in the range of 1 m).
- Massage patch should not contact injury and scar during massage.
- During massage, two massage patches should not be placed around the projection area of the heart at the same time, and it is necessary to make sure that there is no current flowing through the heart no matter how the massage patches are placed.
- It is recommended that when several users use the product, please clean the massage pain relief patch on the product before use.
- When using the product, the massage pain relief patch should fully and evenly contact the skin, otherwise there is a risk of burns.
- The massage patch of this product has a heating function. When the heating function is enabled, the hot compress application on the same part for a long time (more than 20 minutes) may lead to low temperature scalds.
- Discontinue use and consult your doctor if you feel discomfort during or after use.
- This product contains lithium battery. Do not disassemble, modify or soak the product or throw it into fire. Otherwise, the product may be damaged and exploded.
- Do not use the product near the heart, above the neck, head, mouth, private parts or skin disease areas, and do not apply it to the soles of both feet at the same time.
- Do not compress the internal organs during use, otherwise accidents or discomfort may occur.
- Do not use it with other medical devices or multiple sets of this product at the same time, otherwise it may cause nausea or discomfort.
- Do not use it for purposes other than massage as stipulated in the Operation Manual, otherwise accident or malfunction may occur.
- Do not use massage patch with defective or damaged surface, otherwise accident or malfunction may occur.
- Do not disassemble, repair or transform the product without permission. Otherwise, it may lead to failure or electric shock. Read the manual carefully before the use of the product.
- Do not use the product when you carry with electronic instruments, otherwise the electronic instruments may fail (such as: the time of the electronic watch may be disordered).
- If operation failure or other abnormal conditions occur, immediately cut off the power supply and remove the battery. Otherwise, high temperature or short circuit may occur, resulting in an accident. At the same time, contact the customer service personnel for after-sales consultation.

- Do not remove or plug the USB charging cable with wet hands; otherwise, electric shock or accident may occur.
- Do not apply the product to children or people who are unable to express their consciousness, otherwise accidents or discomfort may occur.
- Do not use the product in a place in high humidity, such as bathroom, or use the product while taking a shower, otherwise you may get a strong electric shock and the product may be damaged.
- Do not use it while sleeping, otherwise the product may fail, and the massage patch may be moved to other parts and cause discomfort.
- Do not use it while driving, otherwise it may cause discomfort due to excessive stimulation, resulting in accidents.
- Do not apply the massage patch on the affected area for a long time and use the medium-frequency mode, otherwise it may cause dermatitis.
- During massage, do not touch the massage patch with the metal part of a belt, necklaces and other items, otherwise you may get a strong electric shock.
- When using the massage function, do not apply it in the same part for a long time (more than 30 minutes), otherwise it may cause the muscle fatigue of the application part and physical discomfort.
- Please wear the product correctly: wear it symmetrically on the neck and shoulder, otherwise there may be accidents or physical discomfort.
- The surface of massage patch may become heating, so the people who are not sensitive to heat must be careful when using it.
- After the product works for a long time, the surface will heat up, which is a normal phenomenon.
- Ensure that the instrument button works properly before using it.
- When the product is used for the first time, or the battery power is lost after the product is idle for a long time, please charge the product according to the Operation Manual to ensure that it can work normally.
- Please dispose of this product in accordance with the relevant environmental protection regulations of the city when scrapping it, otherwise it may cause environmental pollution.

Caution

- Read the Operation Manual before using the product.
- If you feel unwell due to the abnormality of this product, please stop using it immediately.
- Consult your doctor and follow the advice of the doctor if you have any discomfort on your body or skin.
- If you want to move the massage patch to another position or other person during the massage, please make sure that the medium frequency is disabled before moving, otherwise you may get a strong electric shock.
- Patients with a history of cerebrovascular accidents should not place the massage patch on the head.

Contraindications:

- Patients with malignant tumors
- People in spasticity after a cerebrovascular accident
- People with high fever (38°C or higher)
- Patients with metal implants
- Patients with spastic paralysis higher)

- Patients with coagulation disorder
- People with anemia
- People with severe cognitive impairment
- People with abnormal neurological stress
- Cardiac pacemaker users
- Eyes
- Patients with infectious diseases
- Near the heart of the human body
- Patients in unconsciousness with contracture deformity of limb bone joint
- Pregnant and parturient women
- Patients with acute suppurative
- Abnormal blood pressure (above 180mmHg)
- Patients with acute [painful] medical conditions, such as spinal fractures, sprains and strains

Please consult your doctor before use it in the following situations, otherwise accidents or discomfort may occur:

- Patients who are under the care of their doctor and feel physically abnormal
- Patients with infectious diseases
- Patients with heart and brain nerve abnormalities
- Patients with impaired skin perception or skin abnormalities
- Patients who are taking drugs
- Patients who cannot freely express themselves
- Patients with perceptual disturbances due to high peripheral circulation disorders such as diabetes
- Patients with abnormal blood pressure

Danger

- Do not use with the following medical electronic devices to avoid injury
- Implanted medical electronic devices such as cardiac pacemakers
- Artificial cardiopulmonary and other life-sustaining medical electronic instruments
- Electrocardiograph and other medical electronic instruments
- Gas stations, high dust workshops, gas stations and other flammable and explosive places

Do not use or store this product in the following places:

- Place with direct sunlight
- Place in high temperature and humidity
- Place with water
- Place with much dust
- Place near fireworks
- Place vulnerable to vibration impact
- Place with strong electromagnetic field