

ACTIVATE YOUR RESULTS

ESSENTIAL FIRST STEP: MICROCURRENT ACTIVATOR

Always prime your skin with NuFACE® Microcurrent Activators before treatment. Each ionized formula is specially formulated to deliver microcurrent to the skin at every level for optimal results and is supercharged with active hydrating or anti-aging ingredients.



FIX® Line Smoothing Serum

Diminishes the look of fine lines and wrinkles, while tightening and brightening for a radiant appearance.

N Ū F A C E®




Follow @MyNuFACE for Aesthetician Tips, Demos, and Live Q+A.
Share your routine or results with #MyNuFACE

75180



N Ū F A C E®
QUICK START





BEFORE YOU POWER UP

CHARGE DEVICE FOR 8 HOURS

Before your first treatment, charge your new NuFACE FIX® Line Smoothing Device and read the enclosed user manual.

SEE HOW EASY TO USE!

Get more out of your new device by visiting [useNuFACE.com](https://www.useNuFACE.com) to watch step-by-step tutorials and learn more.

GETTING STARTED

4 SIMPLE STEPS

1] PREP

Cleanse skin using an oil-free cleanser. This will wash away any excess oils that could otherwise block or prevent an effective microcurrent treatment.

2] ACTIVATE

Apply a mask-like layer of FIX® Serum. Use enough formula to ensure the device glides easily along the skin.

3] SMOOTH

Push the ON/OFF button on your FIX® Device; the white LEDs will illuminate, and the energizing hum will begin. Select a preferred treatment technique. Begin treatment using a quick erasing motion while applying light to medium pressure.

4] FINISH

When complete, massage in remaining FIX® Serum for added skincare benefits.

Visit [useNuFACE.com](https://www.useNuFACE.com) to see how easy to use!

3-MINUTE TREATMENTS



AREAS

TIGHTEN FOREHEAD

SMOOTH FINE LINES
AROUND EYE AREA

REDUCE SMILE LINES

ENHANCE LIP-AREA
FULLNESS

TECHNIQUES

HORIZONTAL
GLIDING MOTION



VERTICAL
GLIDING MOTION



CIRCULAR
MOTION



Optimal results are achieved with regular use.

