

LIFELINKcentral™ AED Program Manager

Contents

- Instructions to log an inspection for one AED at a time
- Instructions to log an inspection for all AEDs at once



To log an inspection for one AED at a time:

1 Visually confirm your AED(s) are “OK” and the battery shows two or more bars

2 Log in

Log into LIFELINKcentral AED Program Manager.

3 Inspections

From “Equipment,” click “Manage Equipment” and then “Inspections” to see the inspection schedules for those AEDs for which you are responsible.

The screenshot displays the LIFELINKcentral AED Program Manager interface. The top navigation bar includes the Stryker logo, the product name, a demo account status, and a user welcome message. A left sidebar contains navigation links: Alerts, My Account, Sites (3), Equipment (10), Manage Equipment, a status summary (AEDs Ready 8 of 10), and a People section (6). The main content area is titled 'Equipment' and features a search bar. A dropdown menu is open under the 'Equipment' header, showing options: Inspections, Device Issues, Custom Equipment, Accessories, and AED Usage Log. The 'Inspections' option is highlighted. Two equipment entries are visible: 'Covington - Remote Office (AA022008)' and 'Redmond - North Building (AA016293)'. The first entry shows a 'LIFEPAK 1000 (IV620357)' with a 'NEEDS ATTENTION (NOTES)' status. The second entry shows a 'LIFEPAK 1000 (IV620359)' with a 'READY' status. Two purple callout boxes with arrows point to the 'Manage Equipment' link in the sidebar and the 'Inspections' option in the dropdown menu.

Manage Equipment

Inspections

4 Inspection Schedules

[Back to contents](#)

Next, click on the inspection schedule to which you would like to navigate. In this example, we will click on "Redmond-North Building."

 **Inspection Schedules**

Enter text to search the list below

Redmond - North Building (AA016293) 

AED Inspection Schedule

Person Responsible	Howard, Ryan
Frequency	Monthly
Last Inspection Date	12/28/2020

Redmond - South Building (AA016294)

AED Inspection Schedule

Person Responsible	Bernard, Andy
Frequency	Monthly
Last Inspection Date	12/28/2020

Navigate to Schedule

5 Select AED

Select the AED for which you are recording the inspection:

LIFEPAK 1000 (IV620357)

Performed By	System
Readiness Time	Not Specified
Location	Medical Clinic

NO DATA AVAILABLE



6 Verify information (Important for first inspection, optional for future inspections)

[Back to contents](#)

Expand each section and verify the information.

Perform Inspection×

LIFEPAK 1000 (IV620357)

EQUIPMENT INFORMATION^

Location

Medical Clinic

Readiness Status

No Data Available

ACCESSORIES∨

INSPECTION CHECKLIST∨

NOTES∨

Submit

Cancel

Expand sections



7 Confirm AED readout

[Back to contents](#)

Proceed to Inspection Checklist and select the matching display readout:

EQUIPMENT INFORMATION

ACCESSORIES

INSPECTION CHECKLIST

Device Present

☒

Visible Damage

☐

Please select the set of indicators exactly as shown on the AED

☒

☐

☐

☐

☐

☐

☐

☐

☐

☐

OK

OK

OK

OK

Readiness Status

Ready

NOTES

Inspection Checklist

- Click the "Submit" button to submit your entries.
- The AED status will turn green for "Ready" confirming a successfully logged inspection:

LIFEPAK 1000 (IV620357)

Performed By

Scott, Michael

Readiness Time

12/28/2020 9:17 AM

Location

Medical Clinic

READY

You will be returned to the list view of the AED Inspection Schedule page. From there, repeat the steps above to perform a maintenance inspection on the next AED in the list.



To log an inspection for all AEDs at once:

[Back to contents](#)

8 Visually confirm the AEDs are “OK” and the battery on each shows two or more bars

9 Log In

Log into LIFELINKcentral AED Program Manager.

10 Inspections

From “Equipment,” click “Manage Equipment” and then “Inspections” to see the inspection schedules for those AEDs for which you are responsible.

Alerts

My Account

Sites3

Equipment10

Manage Equipment

AEDs Ready8 of 10

Ready8

Needs Attention1

No Data Available1

People6

Resource Center

Equipment

Inspect

Enter text to search the list below

Covington - Remote Office (AA022008)

Inspections

Device Issues

Custom Equipment

Accessories

AED Usage Log

LIFEPAK 1000 (IV620357)

Location not specified

CR Plus CHARGE-PAK and 2 Electrodes (1/24/2023)

LIFEPAK 1000 (IV620359)

Security Truck

Battery

Electrodes

LIFEPAK 1000 Battery Non-rechargeable (4/30/2023)

Adult ECG Electrodes (LIFEPAK 1000) (7/26/2023)

Redmond - North Building (AA016293)

LIFEPAK 1000 (IV620357)

Location not specified

NO DATA AVAILABLE

Manage Equipment

Inspections

11 Inspection Schedules

[Back to contents](#)

Next, click on the Inspection Schedule to which you would like to navigate. In this example, we will click on “Redmond-North Building” to log an inspection for all AEDs affiliated with that schedule:



Inspection Schedules

Enter text to search the list below

Redmond - North Building (AA016293)

Navigate to Schedule

AED Inspection Schedule

Person Responsible	Howard, Ryan
Frequency	Monthly
Last Inspection Date	12/28/2020

Redmond - South Building (AA016294)

AED Inspection Schedule

Person Responsible	Bernard, Andy
Frequency	Monthly
Last Inspection Date	12/28/2020



12 AED Inspection Schedule

[Back to contents](#)

Once you are on the landing page for the “Redmond-North Building” Inspection Schedule, click the “Inspect All” button to log a check for all AEDs at once.

The screenshot shows the Stryker LIFELINKcentral AED Program Manager interface. The top navigation bar includes the Stryker logo, LIFELINKcentral AED PROGRAM MANAGER, and a demo customer account for Redmond - North Building. A user profile for Scott, Michael is visible in the top right. The left sidebar contains navigation links for Alerts, My Account, Sites (2), Equipment (4), Manage Equipment, and a summary for AEDs Ready (3 of 4). The main content area is titled 'AED Inspection Schedule' and features an 'Inspect All' button, a search bar, and a 'Complete Inspection' button. Below these are three AED entries: LIFEPAK 1000 (IV620359), LIFEPAK CR Plus (PC12345678), and samaritan PAD 450P (18G00002368). Each entry shows details like 'Performed By', 'Readiness Time', and 'Location'. The first entry has a 'NO DATA AVAILABLE' status, while the others are marked 'READY'. A purple callout box labeled 'Inspect All' points to the 'Inspect All' button in the top right of the main content area.

13 Confirm

The following prompt will appear, asking you to confirm you would like to update the readiness status of all AEDs associated with this schedule to “Ready.” Click “Yes” to complete the check:

The screenshot shows a 'Complete Inspection' dialog box. The text inside reads: 'The readiness status of all Equipment not previously recorded today will be set to 'Ready' with today's date. Do you want to proceed?'. There are two buttons at the bottom: 'Yes' and 'No'. A purple callout box labeled 'Click Yes' points to the 'Yes' button.

14 Completed

[Back to contents](#)

You will be routed to the list view for all AEDs on that Inspection Schedule. The maintenance inspection for all AEDs listed in this schedule is complete, and all AEDs will show green for “Ready” with today’s date; in this example, 12/28/2020 at 9:23 AM:

LIFEPAK 1000 (IV620357)	
Performed By	Scott, Michael
Readiness Time	12/28/2020 9:23 AM
Location	Medical Clinic
READY	
LIFEPAK 1000 (IV620359)	
Performed By	Scott, Michael
Readiness Time	12/28/2020 9:23 AM
Location	Security Truck
READY	



For further information, please contact Stryker at 800 442 1142 (U.S.), 800 668 8323 (Canada) or visit our website at strykeremergencycare.com

Emergency Care

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: HeartSine, LIFELINKcentral, Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker’s trademark or other intellectual property rights concerning that name or logo.

GDR 3336346_C
Copyright © 2021 Stryker

Distributed by:
Physio-Control, Inc.
11811 Willows Road NE
Redmond, WA, 98052 U.S.A.
Toll free 800 442 1142
strykeremergencycare.com

Distributed in Canada by:
Stryker Canada
2 Medicorum Place
Waterdown, Ontario
L8B 1W2
Canada
Toll free 800 668 8323

