
Digital Thermometers

User Guide

Models: 016-638, 016-639

Convenient colour-coded design to support the prevention of cross-contamination.

Features:

- For use orally (blue), rectally (red) or underarm (red recommended)
- Accurate within 0.1°C / 0.2°F
- Gentle-flex tip for added comfort & safety
- Displays in Celsius
- Memory recall / display of last reading
- Beep signals peak temperature
- Waterproof / washable for easy cleaning
- Auto shut-off for extended battery life
- Includes storage case

Thank you for purchasing the Physio Logic® Accuflex® Pro Thermometer. Please read the following instructions for safe use and accurate results.

Your thermometer features a comfortable, soft flexible tip, which conforms to the mouth during temperature readings. It is enhanced with beeper and memory. The following instructions will cover your thermometer’s features in further detail.

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Precautions:

The performance of the device may be impacted should one or more of the following occur:

- Operation outside the manufacturer’s stated temperature and humidity range
- Storage outside the manufacturer’s stated temperature and humidity range
- Mechanical shock (for example, drop test) or degraded sensor
- Patient temperature is below ambient temperature

Cleaning instructions

Before and after each use, clean the thermometer with a soft cloth and 70% isopropyl alcohol diluted with water, or wash with COLD soapy water. An intact unit can be submersed in water or alcohol. However, do not soak in either for an extended period of time.

1. The unit must never be sterilized by boiling, gas, or steam.
2. The use of any other solution or cleaning method beyond the above mentioned may result in reliability problems with the thermometer and void its warranty.

How to take the patient’s temperature

1. Clean and disinfect the thermometer according to the instructions above.
2. Decide whether the measurement will be oral, rectal or underarm.

Note: Measurement time may vary based on the selected site.


Taking an oral temperature

(Blue model recommended)

Note: Drinking hot or cold fluids, exercising, smoking, or performing other activities may raise or lower body temperature. It is important to relax for approximately 5 minutes, mouth closed, prior to taking a reading.

1. Place the probe tip under the tongue as near as possible to the heat pocket, close the patient’s mouth, and avoid moving the unit.
2. Wait a few seconds and then press the ON/OFF button (the unit will beep). The thermometer will first show all displays, then 36.5 °C. The °C symbol will start flashing, which means it is now taking temperature.
3. When the peak temperature has been reached, the °C will stop flashing. The thermometer will beep four times quickly and display the temperature. The display will take a second or two to stabilize and show the final reading. Remove after temperature reading has stabilized.

Note: Temperature reading may increase – even after the 4 beeps – if the thermometer is moved into a warmer heat pocket. If this happens, allow a few seconds for the thermometer to re-stabilize and show final reading. For an accurate reading, the patient’s mouth must remain closed during the measurement. Opening the mouth could result in an extended measurement time and may affect the reading.

Taking a rectal temperature

(Red model recommended)

Note: Rectal temperatures are an appropriate and reliable method for infants and small children. Otherwise, this method should only be used when it is impossible or impractical to take an oral or axillary temperature.

1. Apply a water-soluble lubricant to the tip and insert the probe to a MAXIMUM of 1/2” (1.25 cm) into the rectum.
2. Wait a few seconds and then press the ON/OFF button (the unit will beep). The thermometer will first show all displays, then 36.5 °C. The °C symbol will start flashing, which means it is now taking temperature.
3. When the peak temperature has been reached, the °C will stop flashing. The thermometer will beep four times quickly and display the temperature. The display will take a second or two to stabilize and show the final reading. Remove after temperature reading has stabilized.

Note: Temperature reading may increase – even after the 4 beeps – if the thermometer is moved into a warmer place. If this happens, allow a few seconds for the thermometer to re-stabilize and show final reading.

Note: Rectal temperatures are generally 0.6°C (1.0°F) higher than oral readings.

Taking an axillary temperature (under the arm)

(Red model recommended)

1. Wipe the underarm with a dry towel.
2. Place the probe tip under the arm so the tip is touching the skin, with the thermometer perpendicular to the body.
3. Position the arm across the chest so the probe tip is well covered by the arm. This also ensures that the probe is not affected by the room’s air.
4. Wait a few seconds and then press the ON/OFF button (the unit will beep). The thermometer will first show all displays, then 36.5 °C. The °C symbol will start flashing, which means it is now taking temperature.
5. When the peak temperature has been reached, the °C will stop flashing. The thermometer will beep four times quickly and display the temperature. The display will take a second or two to stabilize and show the final reading. Remove after temperature reading has stabilized.

Note: Axillary temperatures are generally 0.6°C (1.0°F) lower than oral readings.

To turn the unit off

Press the ON/OFF button. However, the unit will automatically shut off after approximately six minutes if it is not turned off manually.

To recall memory

To recall memory of the last reading, hold the ON/OFF button for approximately 2 seconds when turning the thermometer on.


General care

1. Clean the thermometer before and after each use.
2. Store the unit in its protective case when it is not in use.
3. Do not bite the probe.
4. Do not store the unit where it will be exposed to direct sunlight, dust, or humidity. Avoid extreme temperatures.
5. Avoid dropping the thermometer or subjecting it to strong shock.
6. Do not attempt to disassemble the unit, except to replace the battery.
7. The thermometer should only be used under the supervision of an adult.
8. Do not walk, run or talk during temperature taking.

Changing the battery

The battery in your digital thermometer needs to be replaced when the “” symbol appears in the middle right of the display. Replace it with a 1.55V, AG3 alkaline or SR41 silver oxide battery.

1. Slide the battery cap away from the unit.
2. Using a non-metal pointed tool, carefully remove the battery and insert the new battery with the positive (+) side facing the display side of the unit.
3. Carefully slide the battery cap back into place, making sure that the cap fits snugly to protect the battery compartment from moisture.

Note: Please properly dispose of the battery away from small children and heat.

Specifications

Display:	Liquid Crystal Display
Display resolution:	4 Digits (1/10th of a degree increments)
Measurement range:	32°C (89.6°F) to 43.9°C (111°F)
Accuracy:	± 0.1°C (0.2°F) between 35.0°C (95°F) and 39.0°C (102.2°F) ± 0.2°C (0.4°F) otherwise
Average measurement time:	Oral: 5 seconds, rectal: 10 seconds, axillary: 7 seconds
Tone:	Thermometer beeps when activated and deactivated. A series of four beeps signal when peak temperature is reached
Memory:	Displays last measured temperature
Battery:	1.55V, AG3 alkaline or SR41 silver oxide battery
Battery life:	Approximately 90 hours of continuous operation or 18 months of 10 minute use per day
Automatic shut-off:	After approximately 6 minutes
Ambient operating temperature:	0°C (32°F) to 50°C (122°F) (humidity: 10 % ~ 83 %)
Dimensions:	15.2 cm x 2.15 cm x 1.4 cm
Weight :	16.8 grams with battery inside
Specifications are subject to change without notice. For further questions or comments, visit www.amgphysiologic.com	

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