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Instrutions

Aneroid Sphygmomanometers and Cuffs

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Selection of suitable cuff sizes

A. Nylon-velcro cuffs, Disinfectable one piece cuff

Our nylon-velcro cuffs, Disinfectable one piece cuff have a tufted strip on one side and hooks on the other. This ensures that the cuffs can be easily opened and closed quickly and repeatedly. On all the models (with the exception of **sanaphon®**), the cuffs are calibrated, provided with measuring lines. To make sure you have chosen the right cuff size, check to see whether the white index line is in the range between the arrows after the cuff has been put on. If the index line fails to reach this range, the cuff is too small. If it is beyond the range, the cuff is too large. Precise blood pressure readings can be obtained only if the correct cuff size is used.

Nylon-Velcro cuffs (Two Piece Reusable Velcro Cuff). The following cuff sizes are available for all models (except sphgmotensiophone, sanaphon® and exacta®), see cuff table:

Disinfectable one piece cuff (One Piece Reusable Velcro Cuff One Tube/Two Tube). The following cuff sizes are available for all models (except sphgmotensiophone,sanaphon, r-i-san self-measurement and ex-acta), see cuff table:

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Modelles sanaphon® et r-i-san®

5. bigben Round / Square (all versions) and r-i-san, see cuff table:

sanaphon® and r-i-san® (Two Piece Reusable D-Ring Cuff Self-Mesurament): Measure your arm circumference to make sure that it lies within the range indicated on the cuff. The following sizes are available: child, adult, large adult and thigh. These sizes correspond to the circumstances given in the above table.

B. Cotton-velcro cuffs self-measurement (Two Piece Reusable Cuff sphgmotensiophone): Our cotton-velcro cuffs have a tufted strip on one side and hooks on the other. This ensures that the cuffs can be easily opened and closed quickly and repeatedly. Measure your arm circumference to make sure that it lies within the range indicated on the cuff.

The following cuff sizes are available for the models listed below: R1 shock-proof®, minimus® II, minimus® II, bigben Round / Square (all versions) and r-i-san, see cuff table:

C. Cotton hook cuffs (Two Piece Reusable Hook Cuff): On one side of these cuffs, metal bars have been worked into the fabric covering; metal hooks have been riveted to the fabric on the other side. The metal hooks are inserted into the metal bars in the cuff fabric. Measure the circumference of your arm to make sure that it lies within the range indicated on the cuff.

The following cuff sizes are available for the models listed below: R1 shock-proof®, minimus® II, minimus® II, bigben Round / Square (all versions) and r-i-san, see cuff table:

D. Cotton bandage cuffs (Two Piece Reusable Bandage Cuff): There is a bandage strip and a hook on one side of the bandage cuff. To attach the cuff, simply insert the hook into the bandage strip. Measure the circumference of your arm to make sure that it lies within the range indicated on the cuff.

The following cuff sizes are available for the models listed below: R1 shock-proof®, minimus® II, minimus® II, bigben Round / Square (all versions) and r-i-san®, see cuff table:

- Patent cuff used comfortably, the legs not crossed back and arm supported, center of cuff on the upper arm at the level of the right cardiac atrium; the patient should be relaxed as possible and should refrain from talking during the measuring process; approx. 5 minutes should elapse before the first measurement is taken.
- Close the release by turning the deflation screw clockwise [except in the case of **r-i-san®**].
- Putting on the cuff: Put the cuff on in such a way that the lower edge of the cuff is approx. 2 to 3 cm above the crook of the arm level of the right cardiac atrium or approx. 5 cm above the knee joint. Make sure that the marking sign is placed above the artery. The white index strip should be within the marked index range.
- After putting on the cuff, inflate the cuff to approx. 20 mmHg above the expected systolic blood pressure (Diastole = the upper value with the aid of the ball).
- When measuring blood pressure, the position of the operator is normally in front of the patient or at the patient's side.
- It is recommended to use phase V of the Korotkoff sounds (K5) for auscultatory measurement in adults, to use phase IV of the Korotkoff sounds (K4) for auscultatory measurement in children aged 3 to 12 years, and to use phase V of the Korotkoff sounds (K5) for auscultatory measurement in pregnant female patients, except if the Korotkoff sounds can be heard during deflation of the cuff - if this is the case, K4 is to be used.
- Place the chestpiece of the stethoscope, preferably our model anesthesoph, catalog no. 4177-01 - 4177-05, above the artery underneath the cuff.
- For the units intended for home use, no separate stethoscope is needed, since the chestpiece is integrated in the cuff. If a unit intended for home use is used, the membrane of the chestpiece installed in the cuff must be placed above the artery. The cuff is put on by pulling the free end of the cuff through the metal retainer and closing the cuff with the aid of the hook-and-loop fastener.

A. Nylon-velcro cuffs: Close the cuff using the Velcro hook-and-loop fastener.

B. Bandage cuffs: Fasten the bandage cuff by hooking the hook onto the bandage strap.

C. Hook cuffs: In the case of a hook cuff, the metal hook is hooked onto the small metal rods of the cuff cover.

- Patent cuff blood pressure measurement, open the deflation screw by turning it counterclockwise.
- Close the release by turning the deflation screw clockwise [except in the case of **r-i-san®**].
- Putting on the cuff: Put the cuff on in such a way that the lower edge of the cuff is approx. 2 to 3 cm above the crook of the arm level of the right cardiac atrium or approx. 5 cm above the knee joint. Make sure that the marking sign is placed above the artery. The white index strip should be within the marked index range.
- After putting on the cuff, inflate the cuff to approx. 20 mmHg above the expected systolic blood pressure (Diastole = the upper value with the aid of the ball).
- When measuring blood pressure, the position of the operator is normally in front of the patient or at the patient's side.
- It is recommended to use phase V of the Korotkoff sounds (K5) for auscultatory measurement in adults, to use phase IV of the Korotkoff sounds (K4) for auscultatory measurement in children aged 3 to 12 years, and to use phase V of the Korotkoff sounds (K5) for auscultatory measurement in pregnant female patients, except if the Korotkoff sounds can be heard during deflation of the cuff - if this is the case, K4 is to be used.
- Place the chestpiece of the stethoscope, preferably our model anesthesoph, catalog no. 4177-01 - 4177-05, above the artery underneath the cuff.
- For the units intended for home use, no separate stethoscope is needed, since the chestpiece is integrated in the cuff. If a unit intended for home use is used, the membrane of the chestpiece installed in the cuff must be placed above the artery. The cuff is put on by pulling the free end of the cuff through the metal retainer and closing the cuff with the aid of the hook-and-loop fastener.

5.2. mobile stand version / bigben

1. Please assemble the mobile stand first and thereby take note of the installation instruction attached.

2. After assembling the mobile stand, please screw the device (clockwise) onto the mobile stand.

Adjustment:

By opening the retaining screw the desired height can be adjusted. After adjusting the height the mobile stand must be fastened by the retaining screw again.

5.3. ansthesia model / bigben

Remove the wing nut underneath the cuff basket and take off the wall bracket. Affix the wall bracket to the back of 100% rigid plastic bracket so that the wall bracket is placed on the wall. The wall bracket is at the top of the wall bracket engages the edge of the cuff basket and the bottom part fits over the screw that protrudes from the underside of the cuff basket. Now replace and tighten the wing nut on the protruding screw.

again filling up with blood.

- The blood pressure measurement has been completed.
- We would like to point out to you that a unit intended for home use does not replace regular visits to the doctor and that only the doctor can accurately analyze your measured values.

How to care for the aneroid sphygmomanometer

General information

The goal of cleaning and disinfection of medical products is the protection of patients, users and third persons and conserving the value of the medical products. On account of the product design and the use material, no kind of maintenance or repair work can be fixed. The lifetime of the medical products depends on their function and on an appropriate treatment of the devices. Before returning faulty products for repair they must have gone through the described reprocessing process.

1. Manometer and bulb

Manometer and bulb can be cleaned outside with a humid cloth until optical cleanliness is given.

ATTENTION!

Never place the manometer in liquid! This unit is not approved for automated reprocessing and sterilization. These procedures cause irreparable damage!

2. Cuffs

Cotton and Nylon velcro cuff (latex and latex free)

Cleaning:

After removing the bladder, wipe the nylon-velcro covers with a damp cloth. Alternatively, these can be washed with soap and cold water like all the other cuffs. If you decide on the latter course, rinse the cuffs with clear water afterwards and let them air dry. Wipe the bladder and tubes with a damp cloth.

Disinfection:

After removing the bladder, wash the cuff covers in cold water to which disinfectant has been added. After-removing the cuff, dry only disinfectant with approved efficiency and in accordance with the national state data can be used. The bladder and tubes can be wiped with a cotton cloth moistened with ethanol.

Disinfectable one piece cuff

Cleaning:

The cuff can be wiped with a damp cloth. Alternatively, it can be washed with soap and cold water like all the other cuffs. Please rinse the cuff with clear water afterwards. In addition this cuff can be washed up to 60°C in the washing machine, please ensure that no liquid is remaining in the cuff. This can affect the measurement results negatively and it can damage the manometer technology.

Disinfection:

The cuff can be completely inserted into liquid disinfectant. Only disinfectants with approved efficiency and in accordance with the national standards can be used. Before next use, please ensure that no liquid is remaining in the cuff. This can affect the measurement results negatively and it can damage the manometer technology.

IMPORTANT!

Do not iron nylon-velcro cuffs, Disinfectable one piece cuff. Never expose the cuffs to intensive solar radiation! Never touch the cuff covers or bladders with a sharp instrument, since this could cause damage!

Maintenance The product does not require any maintenance.

Test of precision

Remove the tube from the manometer and hold the manometer in a vertical position. When the pointer stands still at 0 on the scale, the instrument has been adjusted properly. If the pointer is below or above 0, the instrument must be recalibrated. Either take it to an authorised **Riester** dealer or send it to us.

Monitoring of instruments

All countries except for Germany:

The respective legal provisions apply for all countries, except for Germany. The reference manometer, which is used for calibration, must be traceable to national and international measurement standards.

Warning: It is not allowed to make changes to the device.

Technical Data:

Round scale

Measurement conditions: 10°C (50°F) to 40°C (104°F) at a relative air humidity of 85 % (non-condensing)

Storage conditions: -20°C (-4°F) to 70°C (158°F) at a relative air humidity of 85 % (non-condensing)

Versions:

- Type of indication: Aneroid, desk, wall, stand and anaesthetic model
- Round scale
- Scale graduation: Increments of 2 mm Hg
- Range shown on scale: 0 to 300 mm Hg
- Measurement range: 0 to 300 mm Hg
- Pointer movement per 1 mm deflation per second: 1 or 2, depending on the particular model
- Tube adapter: Hand-, Tisch-, Wand-, Stand- und Anästhesiemodell
- Pressure generation: Bulb
- Pressure reduction: Air-release valve that can be regulated

Disposal: You can obtain information on disposal from your appropriate local facility or from your local environmental authority.

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Avvertisements / Contre-indications:

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