



CDC-32041-MPK



4 Fr.

1 Lumen

20 cm catheter length

.018 inch dia. spring-wire guide

Pressure Injectable One-Lumen Midline

Contents:

- 1: TaperFree® Catheter††: 4 Fr. (1.45 mm OD) x 20 cm Pressure Injectable Radiopaque Polyurethane with Extension Line Clamp and Contamination Guard

1: GlideThru™ Peel-Away Sheath: 4 Fr. x 2-3/4" (7 cm) Radiopaque over 4 Fr. Dilator

1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer

1: Injection Needle: SafetyGlide™¹ 25 Ga. x 5/8" (1.60 cm) RW

1: Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW

1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)

1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW

1: Syringe: 5 mL Luer-Lock

2: 10 mL Luer-Lock Pre-Filled Saline Syringe

1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep

2: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®² One-Step Solution with Hi-Lite Orange™ Tint

1: Dust Cap: Non-Vented

1: SharpsAway® Disposal Cup

1: SharpsAway® II Locking Disposal Cup

1: Catheter Trimmer

1: Towel: 24" x 36" (61 cm x 91 cm)

1: Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations

1: Safety Scalpel: #11

1: Sterile Procedure Sign

1: Overbed Sign

1: Medication Label: 1% Lidocaine

1: Patient Information Booklet
- 2: Paper Tape Measure

1: Tourniquet

10: Gauze Pad: 2" x 2" (5 cm x 5 cm)

10: Gauze Pad: 4" x 4" (10 cm x 10 cm)

1: Dressing: Tegaderm®³ I.V. 11.5 cm x 8.9 cm

1: Surgical Apparel: Impervious Gown

1: Surgical Apparel: Mask with Eye Shield

1: Dressing: STATLOCK®⁴ PICC Plus Catheter Stabilization Device

1: Surgical Apparel: Bouffant Cap

1: Surgical Apparel: Mask

1: Tape: Steri-Strip®³

¹A trademark of Becton, Dickinson and Company.
²A registered trademark of CareFusion or one of its subsidiaries.
³A registered trademark of 3M Company.
⁴A registered trademark of C. R. Bard, Inc.

***MR Safe
Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).
Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.

 Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)	Pump Flow Rate†† (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Yellow (18 Ga.)	0.37	1880	9324	5

* Priming volumes are approximate and are done without accessories.
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.
†† Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.
** Pressure injection flow rates are determined using media of 11.8 centipoise (cP) viscosity, with 152 cm administration set tubing.

