


SETUP

- Print the sizing strips on LEGAL SIZE paper:
8.5 X 14 in (216 X 356 mm).
- In print settings, set the "Page Sizing & Handling" to "Actual Size" (100 %).
- Use measuring tape to check your printer's scaling.
- Use scissors to cut out the sizing strips that may fit your dog.
- Have your dog lay on its side.




STEP ONE

- Press the finger tab against your dog's leg behind the metacarpal / metatarsal pad's ("heel pad's") furthest point from the front claws, which is not centered behind the heel pad.
- You should be pressing at Clawgs' heel location, which is the leg's skinniest part directly behind the heel pad.

STEP TWO

- Press the sizing strip's loose end against your dog's front pads.
- Make sure your dog's wrist / ankle is bent as far as possible to its standing position.
- The sizing strip's front edge is located at that size's front claws' opening.




SIZE SELECTION

- Select the longest sizing strip that does not contact your dog's front claws' tips while your dog's wrist / ankle is fully bent to its standing position.
- Ideally, the sizing strip will completely cover your dog's front *pad* but will not contact your dog's front *claws' tips*.

Note: Your dog's front paws and back paws may be different sizes.

Note: If you can't decide between two sizes, select the larger size.






dogsap.com

Printable Strips

Clawgs Sizing Guide



IMPORTANT:

- Hold finger tab against the leg behind the heel pad's furthest point from front claws, which is off center.
- Fully bend the wrist/ankle to its 'standing position'

PRINTER SCALING CHECK:

